

Tourni-Cot® makes it safe and easy to achieve hemostasis on digits.

Emergency physicians, hand surgeons, and podiatrists around the world have used the Tourni-Cot® for more than 20 years.

Simple and Easy to Use

Tourni-Cot® saves valuable time over improvised methods to improve workflow.

Quick Hemostasis

Reliable pressure provides a visible field for effective treatment.

Professional Appearance

Sterile and individually packaged, the latex-free Tourni-Cot® communicates that safety comes first.

Minimize Risks

Bright color-coded rings and tag make it easy to select correct size and almost impossible to forget.

Low Cost

Price comparable to sterile surgical gloves. *Data on file



Method	Easy to Use	Immediate hemostasis	Professional appearance	Safe, even pressure	Brightly colored	20 years proven record	Universal for fingers and toes
Tourni-Cot®	✓	✓	✓	✓	✓	✓	✓
Penrose Drain							
Rolled Surgical Glove							
Rubber Band							
Pneumatic Tourniquet			✓	✓			
Plastic Ring	✓	✓	✓	✓			

The Tourni-Cot® Makes a Difference

- Digit injuries commonly have high complication rates
- 38% of Foreign Bodies are missed on first visit (EM Clinics NA, 2007)
- Facilitates the use of tissue and adhesives, simplifying care and reducing visits

“It is critical to wound management” that “the wound is carefully explored in a bloodless field obtained through careful tourniquet application”





(ACEP, Emergency and Primary Care of the Hand)

“Use only CE marked tourniquets which are labeled and/or brightly coloured to enhance visibility. Surgical gloves should not be used for this purpose”


(UK National Patient Safety Agency, 2009)

Indications


- Digit lacerations and avulsions
- Exploration of foreign bodies and deeper
- Nail and nail bed injuries
- Elective procedures

Product Codes & Sizes (product images not to scale)				
Color				
Size	X- Large	Large	Medium	Small
Code	TXL	TCL	TCM	TCS


Using the Tourni-Cot® is Easy for Fingers and Toes



1) Slide proximally to apply and exsanguinate.



2) Even compression is achieved once in place



3) Remove by pulling tag and severing at the notch in the ring.

Place Your Order Today

To learn more about the Tourni-Cot® or to place an order, contact Mar-Med at **1.800.369.3434** or **info@marmedco.com**

MAR-MED

333 Fuller Ave., NE, Grand Rapids, MI 49503 USA

800.369.3434
616.459.5119 Fax

info@marmedco.com
www.marmedco.com